


**Instruction Manual
and Experiment Guide
for the PASCO scientific
Model SE-9640 Franck-Hertz Tube
Model SE-9641 Franck-Hertz Oven
Model SE-9642 Franck-Hertz Control Unit**

012-03498C
2/90

Franck-Hertz Experiment



Copyright © December 1989

\$5.00



PASCO scientific

10101 Foothills Blvd.

• P. O. Box 619011

• Roseville, CA 95661 USA

Table of Contents

| | Page |
|--|------|
| Equipment Return | .ii |
| Introduction | 1 |
| The Franck-Hertz Apparatus | 2 |
| Important Cautions and Tips | .3 |
| The Experiment: Setup, Operation, and Analysis | 4 |

Please— Feel free to duplicate this manual subject to the copyright restrictions below.

COPYRIGHT NOTICE

The PASCO scientific Model SE-9640, **SE-9641**, and **SE-9642** Franck-Hertz Experiment manual is copyrighted and all rights reserved. However, permission is granted to non-profit educational institutions for reproduction of any part of this manual providing the reproductions are used only for their laboratories and are not sold for profit. Reproduction under any other circumstances, without the written consent of PASCO scientific, is prohibited.

LIMITED WARRANTY

PASCO scientific warrants this product to be free from defects in materials and workmanship for a period of one year from the date of shipment to the customer. PASCO will repair or replace, at its option, any part of the product

which is deemed to be defective in material or workmanship. This warranty does not cover damage to the product caused by abuse or improper use. Determination of whether a product failure is the result of a manufacturing defect or improper use by the customer shall be made solely by PASCO scientific. Responsibility for the return of equipment for warranty repair belongs to the customer. Equipment must be properly packed to prevent damage and shipped postage or freight prepaid. (Damage caused by improper packing of the equipment for return shipment will not be covered by the warranty.) Shipping costs for returning the equipment, after repair, will be paid by PASCO scientific.

Equipment Return

Should this product have to be returned to **PASCO** scientific, for whatever reason, notify **PASCO** scientific by letter or phone **BEFORE** returning the product. Upon notification, the return authorization **and** shipping instructions will be promptly issued.

NOTE: NO EQUIPMENT WILL BE ACCEPTED FOR RETURN WITHOUT AN AUTHORIZATION.

When returning equipment for repair, the **units** must be packed properly. Carriers will not accept responsibility for damage caused by improper packing. To be certain the unit

will not be damaged in shipment, observe the following rules:

1. The carton must be strong enough for the item shipped.
2. Make certain there is at least two inches of packing material between **any** point on the apparatus and the inside walls of the carton.
3. Make certain that the packing material *can* not shift in the **box**, or become compressed, thus letting the instrument come in contact with the edge of the **box**.

Introduction

In 1914, J. Franck and G. Hertz performed a landmark experiment, providing important empirical support for Max Planck's quantum theory and for the model of the atom suggested by Niels Bohr. Through the study of collisions between electrons and gas molecules, Franck and Hertz demonstrated that energy is indeed quantized in atomic interactions. For this work, they shared the 1925 Nobel Prize in physics. **Because** of the importance of their experiment, and because of its conceptual simplicity, it has become an important experiment in the undergraduate laboratory.

A simplified **diagram** of the Franck-Hertz experiment is shown in Figure 1. In an oven-heated vacuum tube containing mercury gas, electrons **are** emitted by a heated cathode, and are then accelerated toward a grid that is at a potential V_a relative to the cathode. Just beyond the grid is an anode, which is at a slightly lower potential than that of the grid ($V_p = V_a - \Delta V$; where V_p is the anode potential, and V_a is the grid potential,).

If the accelerated electrons have sufficient energy when they reach the grid, some of them will pass through and will reach the anode. They will be measured **as** current I_c by the ammeter. If the electrons don't have sufficient energy when they reach the grid, they will be slowed by ΔV , and will fall back onto the grid.

Whether the electrons have sufficient energy to reach the anode depends on three factors: the accelerating potential (V_a), the opposing potential (ΔV), and the nature of the collisions between the electrons and the gas molecules in the **tube**. **As** long **as** the electron/molecule collisions are elastic, the collector current depends only **on** V_a and ΔV , because the electrons **lose** no energy to the gas. However, Franck and Hertz discovered that I_c went through a series of maxima and minima **as** the accelerating potential was varied. This implies that the gas molecules absorb energy from the electrons only at specific electron energies (resonant energies).

For example, the first excited state of mercury is 4.9 eV above the ground state. This is the minimum amount of energy that the mercury atoms can absorb in collisions with the accelerated electrons. When V_a is less than 4.9 V, the electron/molecule collisions **are** therefore elastic, so the electrons lose no energy to the gas and arrive at the grid with kinetic energy equal to eV_a . If $V_a > \Delta V$, the electrons have sufficient energy to overcome the opposing voltage, and many of them will pass through the grid and reach the anode, to be measured **as** current. However, if V_a equals 4.9 V, the electrons gain enough kinetic energy to collide inelastically with the mercury atoms just as they reach the grid. In these collisions, the mercury atoms absorb the **full** 4.9 eV




Figure 1 Simplified Diagram of Franck-Hertz Experiment

carried by the electrons. The electrons no longer have sufficient energy to overcome ΔV , and they fall back onto the grid. I_c is then a minimum.

As V_a is **raised** beyond 4.9 V, I_c increases again. However, when V_a reaches 9.8 V, the electrons can lose all their energy to the gas molecules in two collisions. One collision is likely **to** occur midway between the cathode and the grid, the other will occur just as the electrons reach the grid. Again, the electrons lose **all** their kinetic energy in the inelastic collisions, so they fall back onto the grid, and I_c again falls to a minimum. Because of multiple inelastic collisions between the accelerated electrons and the mercury atoms, current minima are found whenever V_a is a multiple of 4.9 V.

Note: The above description is somewhat simplified. Due to contact potentials, the total energy gain of the electrons is not quite equal to eV_a . Therefore, V_a will be somewhat higher than 4.9 V when the first current minimum occurs. However, the contact potentials are a constant in the experiment, **so** at successive current minima V_a will always be a multiple of 4.9 V higher than it was when the first current minimum occurred.

If the accelerating voltage is too high, the electrons can still gain excessive energy before striking Mercury atoms, even if the temperature is correct, and the same problem *can* occur.

How the Tube *Should* Look During the Experiment

When the tube is properly heated and the accelerating voltage is within the proper range, there will be a deep blue cloud between the cathode and the grid (as opposed to the blue-green glow caused by ionization). **You** should **see** dark bands forming near the grid and moving towards the cathode as the accelerating voltage is increased. These bands are the

regions in the tube where the electrons are colliding inelastically with the gas, raising the mercury atoms to their first level of excitation. The blue light is caused by the emission of a photon as the atoms return to the ground state ($E_{\text{photon}} = 4.9 \text{ eV}$; $\lambda_{\text{photon}} = 254 \text{ nm}$).

The Experiment: Setup, Operation, and Analysis

The Franck-Hertz experiment *can* be performed using three different setups:

- Using the Control Unit and an Oscilloscope
- Using the Control Unit and Voltmeters
- Without the Control Unit (using separate power supplies and an ammeter or voltmeter)

The procedure for each setup is described below. Data tables and information on how to analyze the experiment are provided at the end of the manual.

A Using the Control Unit and an Oscilloscope

- Remove the six thumbscrews that secure the front panel of the Franck-Hertz Oven. Replace the front panel with the Franck-Hertz **Tube** as shown in Figure 4.
- Insert a centigrade thermometer (0-200 degrees) into the hole in the top of the oven. Insert it **so** the tip of the thermometer is near the center of the **tube**.
- Check the label and plug the oven into an appropriate AC power outlet, then turn the thermostat dial to 180 °C. Let the oven heat up for 10-15 minutes. Keep an eye on the thermometer. Do **not let the oven temperature exceed 205 °C**.

- Connect the tube, control unit, and oscilloscope as shown in Figure 4. (The voltmeters that are shown in the figure are not needed when an oscilloscope is used.)
- Switch on the control unit, but only AFTER the tube has warmed up for 10-15 minutes. Set the controls as follows:

| | |
|------------|--|
| Heater: | midrange (about 5.5 V) |
| Amplifier: | Gain—midrange Zero—midrange |
| Va: | Adjust—not applicable Switch—ramp (\rightarrow) |

- Wait about 90 seconds for the cathode to warm up, then adjust the gain controls on the oscilloscope to get a trace as shown in Figure 4. (The x-channel gain should be approximately 0.5 volts/cm; the y-channel gain should be approximately 0.5 volts/cm.)
- From the oscilloscope trace, measure the values of V_a for which the collector current (I_c) is a minimum.
Note: The input to the y-channel of the oscilloscope is proportional to I_c . The input to the x-channel of the oscilloscope is equal to $V_a/10$.
- See Analyzing the Experiment at the end of the manual for details about manipulating and analyzing your data.




Figure 4 Franck-Hertz Experiment Using the Control Unit

B. Using the Control Unit and Voltmeters

1. Remove the six thumbscrews that secure the front panel of the Franck-Hertz Oven. Replace the front panel with the Franck-Hertz **Tube** as shown in Figure 4.
2. Insert a centigrade thermometer (0-200 degrees) into the hole in the top of the oven. Insert it so the tip of the thermometer is near the center of the **tube**.
3. Check the label and plug the oven into an appropriate AC power outlet, then **turn** the thermostat dial to 180 °C. Let the oven heat up for 10-15 minutes. **Keep an eye on the thermometer. Do not let the oven temperature exceed 205 °C.**
4. Connect the **tube**, control **unit**, and voltmeters as shown in Figure 4 (the oscilloscope connections that are shown in the figure **are** not necessary when using voltmeters).
5. Switch on the control unit, but only **AFTER** the **tube** has been allowed to warm up for 10-15 minutes. Initially, set the controls **as** follows:

| | |
|------------|---|
| Heater: | midrange (about 5.5 V) |
| Amplifier: | Gain - maximum gain (fully clockwise) Zero - midrange |
| V_a : | Switch—Man. (→) Adjust - zero (fully counterclockwise) |

6. Allow at least 90 seconds for the cathode to heat up.
7. Adjust the Amplifier Zero so that the amplifier output is zero volts (if it won't go to zero, just adjust it **as** close to zero **as** possible).
8. Increase the accelerating voltage in increments of about 0.5 volts, by turning the V_a Adjust knob clockwise. At each value of V_a , record V_a and the corresponding value for the amplifier output.

Note: The amplifier output is not equal to I_c . It is a voltage that is proportional to I_c .

9. **See Analyzing** the Experiment at the end of the manual for details about manipulating and analyzing your **data**.

C. Without the Control Unit

1. Remove the six thumbscrews that secure the front panel of the Franck-Hertz Oven. Replace the front panel with the Franck-Hertz Tube as shown in Figure 4.
2. Insert a centigrade thermometer (0-200 degrees) into the hole in the top of the oven. Insert it so the tip of the thermometer is near the center of the tube.
3. Check the label and plug the oven into an appropriate AC power outlet, then turn the thermostat dial to 180°C. Let the oven heat up for 10-15 minutes. Keep an eye on the thermometer. **Do not let the oven temperature exceed 205°C.**
4. Connect the tube, power supplies, and ammeter* as shown in Figure 5.

Important: All power supplies must be floating with respect to ground.

6. Set all the power supplies to zero volts. Then turn them on and set the heater voltage to 5.5 volts, the accelerating voltage to zero volts, and the opposing voltage to 1.5 volts. **DO NOT LET THE HEATER VOLTAGE EXCEED 7 VOLTS**, or the tube may be damaged.

17. Increase the accelerating voltage, V_a , in increments of about 0.5 volts. At each value of V_a : record V_a and the corresponding value for the anode current, I_a . Repeat your measurements, but this time **try** to locate the exact values of V_a for which current minima occur. Record these values of V_a and the corresponding values of I_a .
18. **See Analyzing the Experiment** at the end of the manual for details about manipulating and analyzing your **data**.

***Note:** If you don't have a sufficiently sensitive ammeter to accurately measure the anode current, you *can* use a voltmeter instead. The voltmeter must have 0.1 mV resolution and a 10M Ω input impedance. Connect it just as the ammeter is connected in Figure 5. The voltage reading of the voltmeter will then be **equal** to the anode current, in units of 10⁻¹¹ amperes.




Figure 5 Franck-Hertz Experiment Without the Control Unit